



Basic and Applied Scientific Research Centre Dammam 31113, Saudi Arabia		 جامعة الإمام عبد الرحمن بن فيصل IMAM ABDULRAHMAN BIN FAISAL UNIVERSITY مركز البحوث العلمية الأساسية و التطبيقية Basic and Applied Scientific Research Center	
Procedure Title: QUALITY MANUAL ISO/IEC 17025:2017		Document No. Samples -IAU-BASRC-010201-F-06	
Department/Section	Quality Assurance / Lab	Revision Date: 1.1.2022	Next Revision Date: 1.1.2026
		Revision No. 1	Page 1 of 2

REQUEST FOR TGA/ DTA SERVICE

TYPE OF SERVICE:	<input type="checkbox"/> Internal	<input type="checkbox"/> External
-------------------------	-----------------------------------	-----------------------------------

REQUESTER INFORMATION	
Client Name:	
Organization Type: <input type="checkbox"/> College <input type="checkbox"/> Research institute <input type="checkbox"/> University	
Name of organization:	
Department:	
E-mail:	
Mobile:	
Sig.:	Date:

Basic and Applied Scientific Research Centre Dammam 31113, Saudi Arabia		 جامعة الإمام عبد الرحمن بن فيصل IMAM ABDULRAHMAN BIN FAISAL UNIVERSITY مركز البحوث العلمية الأساسية و التطبيقية Basic and Applied Scientific Research Center	
Procedure Title: QUALITY MANUAL ISO/IEC 17025:2017		Document No. Samples -IAU-BASRC-010201-F-06	
Department/Section	Quality Assurance / Lab	Revision Date: 1.1.2022	Next Revision Date: 1.1.2026
		Revision No. 1	Page 2 of 2

SAMPLE INFORMATION
Number of samples:
Required temperature (Celsius):
Sample analysis: <input type="checkbox"/> With Nitrogen gas <input type="checkbox"/> Without Nitrogen gas
Safety precautions: <input type="checkbox"/> Toxic <input type="checkbox"/> Flammable <input type="checkbox"/> Biohazardous <input type="checkbox"/> None
INSTRUCTIONS
1- If the temperature is higher than 600, it must be ensured that the sample does not interact with platinum .
2- If its 600 or less, it must be ensured that the sample does not interact with Aluminum .
3-The results will be sent by e-mail <u>within 10 working days</u>
Would you like to retrieve your samples after completing the analysis: <input type="checkbox"/> Yes <input type="checkbox"/> No
Please collect your samples no later than one week from receiving this e-mail. After a week, the samples will be disposed
Storage requirements for the sample :.....
Acknowledgment: All information filled in by the Client is correct

Samples Supervisor Name:

.....

Sig.:

Date:

